**NGS**

# ***Korea Medical Group***

LabGenomics Co., Ltd

2019

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Disease	Gene	Risk Allele	Ref. Allele	Your Allele	Result	Relative Risk
Liver Cancer	STAT4	G	T	G G	High Risk	2.01
	HLA-DRB1 - LOC107986589	A	G	A A		
	MTCO3P1 - LOC102725019	A	G	A G		
Pancreatic cancer	LOC105370243	G	T	G G	High Risk	1.76
	LINC01394 - LOC105374880	C	T	C C		
Prostate cancer	CASC8, CCAT2	G	T	G G	High Risk	1.9
	RFX6	T	C	T T		
Testis cancer	BAK1	G	A	A G	High Risk	1.56
	KITLG	G	A	G G		
Colorectal cancer	C5orf66	A	C	C A	Average Risk	1.18
	MYRF	G	T	T T		
	CASC8, CCAT2	G	T	G G		
Renal cell carcinoma	SCARB1	T	C	C T	Average Risk	1.19
	LOC105369705	G	A	G G		
Thyroid cancer	PTCSC3 - RN7SKP21	T	C	C C	Low Risk	0.76
	DIRC3	T	C	T T		
Gastric cancer	PSCA	T	C	C C	Low Risk	0.47
	PRKAA1	C	T	T T		
Bladder cancer	CASC11	T	G	G G	Low Risk	0.65

Sample ID 2018-04-22 00013

의뢰기관명 Institution	검사대상자명 Name of Patient	차트번호 Chart No.
Sample	Sample	Sample

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Disease	Gene	Risk Allele	Ref. Allele	Your Allele	Result	Relative Risk
Bladder Cancer	C20orf187	A	C	C C	Low Risk	0.65
	HECTD4	A	G	G G	Low Risk	
Esophageal Cancer	ALS2CR12	A	G	G G	Low Risk	0.22
	TERT	G	A	A A	Low Risk	
Lung Cancer	BPTF	A	G	G A	Low Risk	0.57
	CDKN2B-AS1	T	C	T T	High Risk	
Cerebral aneurysm	RP1	A	G	A A	High Risk	1.43
	LOC101929163, BTNL2	T	C	C T	High Risk	
Osteoarthritis	ALDH1A2	C	G	G C	High Risk	1.45
	PRDM16, LOC105378606	C	T	T C	High Risk	
Migraine	MEF2D	C	T	T T	High Risk	2.06
	FTO	A	T	T T	Average Risk	
Obesity	MC4R	C	T	T T	Average Risk	0.84
	CLOCK	A	G	A A	Average Risk	
Hypertension	UMOD	A	G	A A	Average Risk	0.81
	FGF5	C	T	T T	Average Risk	
Stroke	PITX2 - MIR297	G	A	G A	Average Risk	0.99
	SPSB4	G	A	G G	Average Risk	
Atrial fibrillation	ZFX3	T	C	C T	Average Risk	1.22
	PITX2 - MIR297	C	T	T C	Average Risk	
Chronic obstructive pulmonary disease	KRT18P51 - HHIP-AS1	T	C	T T	Average Risk	0.77
	FAM13A, LOC105377327	C	T	T T	Average Risk	
Asthma	HLA-DQB1	T	C	T C	Average Risk	1.03
	IKZF4	G	T	T G	Average Risk	
Osteoporosis	WHSC1L2P - SOST	T	C	T T	Average Risk	0.76
	MCOM	A	G	G G	Average Risk	

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Registration No. 2018-04-22 99913

의뢰기관명 Institution	검사대상자명 Name of Patient	차트번호 Chart No.
Sample	Sample	Sample

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
Disease	Gene	Risk Allele	Ref. Allele	Your Allele	Result	Relative Risk
Alzheimer's disease	CLU	T	C	C T	Average Risk	1.08
	SORL1	T	C	T T		
Depressive disorder	CNTN5	C	T	C T	Average Risk	0.92
	KSR2	A	G	G G		
Diabetes mellitus type 2	LOC105375716, SLC30A8	C	T	C T	Low Risk	0.67
	KCNQ1	C	T	T T		
Myocardial infarction	AP3D1 - DOT1L	C	A	C A	Low Risk	0.81
	PLCL2	G	A	A G		
Parkinson's disease	MCCC1	G	A	G A	Low Risk	0.64
	HLA-DRA	G	A	A A		

Registration No. 2018-04-22 99913

의뢰기관명 Institution	검사대상자명 Name of Patient	차트번호 Chart No.
Sample	Sample	Sample

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DOB: 1975-06-05

Age/Gender: 43 / F

Patient Information	Number of Fetus	Single	Ultrasonographic Finding	
	Ultrasound Gestational Age	10W 1D	Prenatal Biochemical Screening findings	
	Height/Weight	165 cm / 65 kg		

#### Quality Control / Test Result

Cell free DNA quality	Sequencing quality	Fetal fraction	Standard Material Test result
Pass	Pass	Pass(15.28%)	Pass

#### INTERPRETATION

Low Risk - The chance of the baby having a chromosomal abnormality is very low.

#### RESULT DETAILS

Chromosome Aneuploidies		
Items (Disease Type)	Result	Risk Score
Trisomy 21 (Down Syndrome)	Low Risk	1/12435
Trisomy 18 (Edward Syndrome)	Low Risk	1/38045
Trisomy 13 (Patau Syndrome)	Low Risk	1/12737
Sex aneuploidy (mX, XXY, XXX)	Low Risk	1/5557, Male

**Test Method**

cfDNA is isolated from maternal blood and sequenced using Next Generation Sequencing (NGS) technology. Sequencing data is analyzed using in-house bioinformatics pipeline to identify fetal aneuploidy in the tested chromosomes.

\*This test provides a result only when a sample meets the quality threshold.

**Test Purpose & Limitation**

\*This test is a screening test for T21, T18, T13 and sex chromosome aneuploidy under the consent of the mother.

\*No analysis for sex chromosome aneuploidy is to be performed in case of twin.

\*This test is highly accurate, but not diagnostic. Therefore, if a confirmatory test is required according to the test result or the clinical situation of the mother, an amniocentesis or CVS should be performed.

MomGuard Performance Table		
	Sensitivity	Specificity
Chromosome 21	98.65%	99.94%
Chromosome 18	100%	99.98%
Chromosome 13	100%	99.97%

Overall PPV = 0.8932, NPV=0.9999  
Overall False Positive (FP) % = 0.102  
Overall False Negative (FN) % = 0.009

LabGenomics Co.,Ltd.

Genetic Testing Laboratory: No. 23  
Genetic Research Laboratory: No. 7

Inspector: MJ OH Ph.D.

Analysis Officer: HJ HU Ph.D.

Laboratory Director: SY KIM M.D., Ph.D.

Laboratory Officer: DY CHO Ph.D.

LGT\_1301\_11E1 Rev(2) 18.07.01

### Test Report

Registration No.	E123456789
Collection Date	2018 Y 05 M 13 D
Receipt Date	2018 Y 05 M 20 D
Analysis Date	2018 Y 05 M 22 D
Report Date	2018 Y 05 M 23 D

Department: /

Doctor:

Sample ID:

Address:

QR code

Organization: Sample

Patient: Sample

Chart No.:

DOB: 1977.06.25

AGE/Gender: 40 / F

Patient Information	Number of Fetus	SINGLE	Ultrasonographic Finding	
	Ultrasound Gestational Age	15 W 4 D	Prenatal Biochemical Screening findings	
	Height/Weight	163 cm / 83 kg		

#### Quality Control / Test Result

Cell free DNA quality	Sequencing quality	Fetal fraction	Standard Material Test result
Pass	Pass	Pass (12.06 %)	Pass

#### INTERPRETATION

**Trisomy 21 ( Down Syndrome ) : High Risk**  
Increased risk of Trisomy 21. Genetic counseling and confirmatory testing is required.

Chromosome Aneuploidies			Chromosome Microdeletions	
Tested condition	Result	Risk Score	Tested condition	Result
Trisomy 21 (Down Syndrome)	High Risk	> 9/10	Cri-du-chat	Low Risk
Trisomy 18 (Edward Syndrome)	Low Risk	<1/12,300	1p36 deletion	Low Risk
Trisomy 13 (Patau Syndrome)	Low Risk	<1/11,435	DiGeorge[22q11.2 deletion]	Low Risk
Trisomy 9	Low Risk	<1/10,000	Jacobsen	Low Risk
Trisomy 16	Low Risk	<1/10,000	Prader-Willi/Angelman	Low Risk
Trisomy 22	Low Risk	<1/10,000		
Sex aneuploidy (mX, XXY, XXX)	Low Risk	<1/10,000		

**Test Method**

cfDNA is isolated from maternal blood and sequenced using Next Generation Sequencing (NGS) technology. Sequencing data is analyzed using in-house bioinformatics pipeline to identify fetal aneuploidy in the tested chromosomes.

\*This test provides a result only when a sample meets the quality threshold.

**Test Purpose & Limitation**

\*This test is a screening test for T21, T18, and T13 with the following options under the consent of the mother: additional chromosomal aneuploidies (T9, T16, T22), sex chromosome aneuploidy and selected microdeletion syndromes.

\*This test is highly accurate, but not diagnostic. Therefore, if a confirmatory test is required according to the test result or the clinical situation of the mother, an amniocentesis or CVS should be performed.

MomGuard Performance Table		
	Sensitivity	Specificity
Chromosome 21	98.65%	99.94%
Chromosome 18	100%	99.98%
Chromosome 13	100%	99.97%

Overall PPV = 0.8932, NPV=0.9999  
Overall False Positive (FP) % = 0.103  
Overall False Negative (FN) % = 0.009

LabGenomics Co.,Ltd.

Genetic Testing Laboratory: No. 23-2  
Genetic Research Laboratory: No. 7-2

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QR code

### EnfantGuard2.0 Test Report

Sample ID:	99990	Organization:	SNU hospital
Clinician:	Erin K	Name of Patient:	Amy Lee
Clinical Opinion:	N/A	DOB/Gender:	2017. 08. 08 / F
Department:	/	Chart No.:	13526
Address:		Sample Type:	Capillary tube

**Test Information**


This test is a screening test for chromosome abnormality with neonatal developmental disorder. This test is performed using NGS to analyze more than 6 million sequences.  
This test is only for personal reference but not for clinical diagnosis. For abnormal result subject, appropriate genetic counseling and a confirmation test are required.

**Quality Control**

DNA quality	Sequencing quality	Analysis Quality Control Result
Pass	Pass	Pass


**Test Result**

<b>Test Result</b>	High Risk
<b>Description</b>	chr17:14125001-15775000, 1.65 Mb, deletion Micro-deletion of chromosome 17 was found, related with 17p12 deficiency syndrome.
<b>Comment</b>	<ul style="list-style-type: none"> <li>- This test is a screening test for chromosome abnormality with neonatal developmental disorder.</li> <li>- This test cannot be used to detect other chromosomal defects(balanced translocation, inversion, point mutation, low level mosaicism, etc.).</li> <li>- Possibility of other genetic alterations that are not included in this test can not beruled out.</li> <li>- This test provides the results only if it satisfies the specified qualitystandards.</li> </ul>

Genetic Testing Laboratory : No.23-2  
Genetic Research Laboratory : No.7-2 MU Oh Ph.D.

Inspector :  Oh Ph.D.

Analysis Officer :  SJ No Ph.D.

Laboratory Officer :  DY Cho Ph.D.

Sampling Date	2016 Y 12 M 25 D
Request Date	2016 Y 12 M 25 D
Report Date	2017 Y 01 M 03 D

Name of Organization	Name of Patient	Sample ID
SNU hospital	Jenny P	99990


Chromosome 1

Result

Negative

Chromosome 2

Result

Positive

ChrX:155270560-155270560, 155.27Mb, deletion

Deafness, Dystonia, a and Cerebral Palsy/epilepsy, 17p13.3. Centromeric duplication syndrome, Parkinson type of Early-Onset Parkinson Disease

Chromosome 3

Result

Negative

Chromosome 4

Result

Negative


Chromosome 5

Result

Negative

Chromosome 6

Result

Positive

ChrX:1-155270560, 155.27Mb, deletion

A deletion in chromosome 6 has been found

B-6F, 700, Daewangpangyo-ro, Bundang-gu, Seongnam-si, Gyeonggi-do, Korea 13489

Genetic Testing Laboratory : No.23-2  
Genetic Research Laboratory : No.7-2

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1/3 page

### BRCA1/BRCA2 GENETIC TEST REPORT

Collection Date	
Receipt Date	
Report Date	

Organization		SPECIMEN	
Name	SAMPLE	Specimen #	
Address		Specimen Type	

PHYSICIAN		PATIENT	
Name	SAMPLE	Name	
Contact Info		Date of Birth /Gender	
Clinical comment		Medical Chart #	

#### SUMMARY RESULT

[A p|P]athogenic variant[|s] [was|were] detected in the [BRCA2 gene|BRCA1 gene|BRCA1 and BRCA2 genes]

[VUS]  
X variant[s] of uncertain significance (VUS) [was|were] detected.

[Comment]  
This variant[s] [is|are] clinically significant and [is|are] associated with an increased cancer risk. In this case, **genetic counseling and a BRCA1,2 genetic test of your family members are recommend**. This result does not mean that you have a diagnosis of cancer or that you will definitely develop cancer in your lifetime. Your actual risk may be different based on other genetic and non-genetic factors.

2/3 page

### BRCA1/BRCA2 GENETIC TEST REPORT

Collection Date	
Receipt Date	
Report Date	

#### DETAILED RESULT

##### [Pathogenic variant]

Gene	Classification	Zygotity	Variant Detected	Amino Acid Change	dbSNP
BRCA2	Pathogenic	Hetero	c.5074G>C	p.Asp1692His	rs80187739

A pathogenic variant ([homozygous | heterozygous] c.XXXXXX (p.XxxXXXXxx)) was detected in the [exon X| intron] of the BRCA2 gene ([NM\_000059|NM\_007294]).

This [nonsense | missense | silent | insertion | deletion | frameshift] mutation has been previously reported [once|twice|X times] in our database before (NM\_00XXXX). Variants detected within each gene are reported and classified according to our internal criteria with reference to ACMG guidelines

##### [VUS]

Gene	Classification	Zygotity	Variant Detected	Amino Acid Change	dbSNP
BRCA2	VUS	Hetero	c.29C>A	p.Thr10Lys	rs1057519494


X VUS (BRCA2 c.XXXXXX (p.XxxXXXXxx); BRCA2 c.XXXXXX (p.XxxXXXXxx)) with unknown relationship with disease [was|were] detected

In addition, X benign variant[s] which [has|have] no impact on health [was|were] detected in the [BRCA1 gene|BRCA2 gene|BRCA1,2 genes].

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검체정보: FFFE (Specimen Information)  
 조직정보: 대장암 (Tissue Information)  
 Tumor purity: 80% (Estimated by Pathologist)  
 병리진단소견: Adenocarcinoma (Diagnosis)

의뢰기관명: FFFE (Institution)  
 진료과/병동: (Department)  
 의뢰의사명: (Name of Doctor)  
 차트번호: (Chart No.)  
 검사대상자명: (Name of Patient)  
 생년월일: (Date of Birth)  
 나이/성별: (Age/Gender)

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**Negative** (Tumour variants not found<sup>4</sup>)

Gene	Variant	Missed Depth CDS
BRCA2	Wild Type	1589x
BRAF	Wild Type	888x
EGFR	Wild Type	773x
ERBB2	Wild Type	613x
NRAS	Wild Type	1012x

**CancerSCAN Variants: 29**  
 - Annotated Variants: 3  
 - Known Variants: 2  
 - Novel Variants: 24

Estimated tumor purity (Bioinformatics): 73%  
<sup>1</sup> For the variant classification, please see the supplementary report.  
<sup>2</sup> The variant frequency observed from TOGA datasets, 7,225 clinical samples of 24 tumor types (Ref: [www.cbioportal.org](http://www.cbioportal.org)).  
<sup>3</sup> The variant frequency observed from MSK-IMPACT datasets, 10,945 clinical samples of 20 tumor types (Ref: [ref.mskcc.org](http://ref.mskcc.org)).  
<sup>4</sup> In the genes (i.e. TP53, BRAF, KRAS, EGFR, BRCA1, BRCA2 and ERBB2) are fully covered (mean depth > 200x) and tumor hotspot positions were not found in these genes, these are recorded as tumor-mutation free.  
<sup>5</sup> Known and novel variants are listed in supplementary report.

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**상세결과 I | Result Detail - Part A**

**A.Variant Annotation (Table view)**

Variant	Disease	Class	CIVIC Variant	Drug	Clinical significance
KRAS G12D	Colorectal Cancer	1	EXON 2 MUTATION	Cetuximab	Resistant
			MUTATION	Cetuximab	Resistant
	Colorectal Cancer	2	MUTATION	Platumumab	Resistant
			EXON 2 MUTATION	Bevacizumab	Resistant
	Colorectal Cancer	2	EXON 2 MUTATION	Cetuximab	Resistant
			G12G13	Cetuximab	Resistant
			MUTATION	Cetuximab	Sensitive
			MUTATION	Chemotherapy	Resistant
			EXON 2 MUTATION	EGFR Inhibitor	Resistant
			EXON 2 MUTATION	Immunecan	Sensitive
Non-small Cell Lung Carcinoma	3	EXON 2 MUTATION	Regorafenib	Sensitive	
		MUTATION	Regorafenib	Sensitive	
		MUTATION	RO4987655	no info	
		EXON 2 MUTATION	Selumetinib (AZD6244)	Sensitive	
Non-small Cell Lung Carcinoma	3	MUTATION	Selumetinib (AZD6244)	Sensitive	
		G12	(Prognostic)	no info	
		G12D	(Prognostic)	Better Outcome	
Acute Leukemia	4	G12	(Diagnostic)	Positive	
Endometrial Cancer	4	MUTATION	Ridorsolimus	Sensitive	
Hepatocellular Carcinoma	4	MUTATION	Temsirolimus	Sensitive	
		MUTATION	Regorafenib	Sensitive	
		MUTATION	Sorafenib	Sensitive	

<sup>4</sup> Annotation information of variants was from CIVIC (<https://civic.genome.wustl.edu>, version 2017-09-01)

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**B.Variant Annotation (Description)**

**EGFR L858R**

EGFR L858R has long been recognized as a functionally significant mutation in cancer, and is one of the most prevalent single mutations in lung cancer. Best described in non-small cell lung cancer (NSCLC), the mutation seems to confer sensitivity to first and second generation TKIs like gefitinib and neratinib. NSCLC patients with this mutation treated with TKIs show increased overall and progression-free survival, as compared to chemotherapy alone. Third generation TKIs are currently in clinical trials that specifically focus on mutant forms of EGFR, a few of which have shown efficacy in treating patients that failed to respond to earlier-generation TKI therapies.

**Class 1 : Non-small Cell Lung Carcinoma - L858R - Erlotinib - Sensitive**  
 Ref: Khoo et al., 2014, Clin Oncol. PMID:24680989  
 U.S. Food and Drug Administration approval summary: Erlotinib for the first-line treatment of metastatic non-small cell lung cancer with epidermal growth factor receptor exon 19 deletions or exon 21 (L858R) substitution mutations.

On May 14, 2013, the U.S. Food and Drug Administration approved erlotinib (Tarceva) for the first-line treatment of patients with metastatic non-small cell lung cancer (NSCLC) whose tumors have epidermal growth factor receptor (EGFR) exon 19 deletions or exon 21 (L858R) substitution mutations.

**Class 2 : Lung Adenocarcinoma - L858R - Afatinib - Sensitive**  
 Ref: Sequist et al., 2013, J. Clin. Oncol. PMID:23819960  
 Phase III study of afatinib or cisplatin plus pemetrexid in patients with metastatic lung adenocarcinoma with EGFR mutations.

A phase III clinical trial (NCT00949650) found that median progression-free survival among patients with exon 19 deletion or L858R EGFR mutations (n = 308) was 13.6 months for afatinib and 6.3 months for chemotherapy (HR, 0.47; 95% CI, 0.34 to 0.65; P < 0.001).

Ref: Yang et al., 2012, Lancet Oncol. PMID:22403297  
 Afatinib for patients with lung adenocarcinoma and epidermal growth factor receptor mutations (LUX-Lung 2): a phase 2 trial.

In a phase 2 study of patients with lung adenocarcinoma (stage IIB with pleural effusion or stage IV) and EGFR mutations, treated with afatinib were assessed by objective response. 129 patients were treated with afatinib. 66% of the 106 patients with less common activating EGFR mutations (deletion 19 or L858R) had an objective response compared to 39% of 23 patients with less common mutations.

Ref: Wu et al., 2014, Lancet Oncol. PMID:24439029  
 Afatinib versus cisplatin plus gemcitabine for first-line treatment of Asian patients with advanced non-small-cell lung cancer harbouring EGFR mutations (LUX-Lung 6): an open-label, randomised phase 3 trial.

Afatinib is an irreversible covalent inhibitor of EGFR (second generation). This Phase III clinical trial (LUX-Lung 6; NCT01121383) was performed in Asian patients with EGFR mutant advanced NSCLC. 364 eligible patients with EGFR mutations were assigned to afatinib (n=242) or gemcitabine and cisplatin (n=122) treatment. The trial observed significantly longer median progression-free survival with afatinib vs. gemcitabine and cisplatin treatment (11.0 vs. 5.6 months). Afatinib/Chemotherapy group composition: 51.250.8% del 19; 38.97.7% L858R; 10.811.5% Unknown.

**Class 2 : Non-small Cell Lung Carcinoma - L858R - Gefitinib - Sensitive**  
<sup>4</sup> Annotation information of variants was from CIVIC (<https://civic.genome.wustl.edu>, version 2017-09-01)

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#### Negative (Tumour variants not found<sup>4</sup>)

Gene	Variant	Mean.Depth.CDS
ALK	Wild Type	982x
BRAF	Wild Type	531x
BRCA1	Wild Type	1,148x
BRCA2	Wild Type	959x
ERBB2	Wild Type	815x
IDH1	Wild Type	997x
IDH2	Wild Type	1,054x
KIT	Wild Type	1,054x
MYC	Wild Type	1,082x
MYCN	Wild Type	908x
NRAS	Wild Type	1,204x
PDGFRA	Wild Type	485x

Estimated tumor purity (Bioinformatics): 13%

**CancerSCAN Variants : 21**  
**- Annotated Variants : 3**  
**- Known Variants<sup>5</sup> : 9**  
**- Novel Variants<sup>5</sup> : 9**

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**Negative** (Tumour variants not found<sup>4</sup>)

Gene	Variant	Mean.Depth.CDS
BRCA2	Wild Type	1589x
BRAF	Wild Type	555x
EGFR	Wild Type	772x
ERBB2	Wild Type	812x
NRAS	Wild Type	1052x

**CancerSCAN Variants : 29**  
**- Annotated Variants : 3**  
**- Known Variants : 24**

Estimated tumor purity (Bioinformatics): 78%

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
Variant	Disease	Class	CIViC Variant	Drug	Clinical significance
KRAS/G12D	Colorectal Cancer	1	EXON 2 MUTATION	Cetuximab	Resistant
			MUTATION	Cetuximab	Resistant
			MUTATION	Panitumumab	Resistant
Colorectal Cancer	Colorectal Cancer	2	MUTATION	Bevacizumab	Resistant
			EXON 2 MUTATION	Cetuximab	Resistant
			G12/G13	Cetuximab	Resistant
			MUTATION	Chemotherapy	Resistant
			EXON 2 MUTATION	EGFR Inhibitor	Resistant
			EXON 2 MUTATION	Irinotecan	Sensitive
			EXON 2 MUTATION	Selumetinib (AZD6244)	Sensitive
			G12D	(Prognostic)	Better Outcome
Non-small Cell Lung Carcinoma	Non-small Cell Lung Carcinoma	3	MUTATION	Docetaxel	Resistant
			MUTATION	Selumetinib (AZD6244)	Resistant
Acute Leukemia	Acute Leukemia	4	G12	(Diagnostic)	Positive
Hepatocellular Carcinoma	Hepatocellular Carcinoma	4	MUTATION	Refametinib	Sensitive
			MUTATION	Sorafenib	Sensitive
Lung Adenocarcinoma	Lung Adenocarcinoma	4	MUTATION	(Prognostic)	Poor Outcome
Lung Cancer	Lung Cancer	4	G12D	Gefitinib	Resistant
			G12D	(Diagnostic)	Positive

상세결과 | Result Detail - Part A

A. Variant Annotation (Table view)

Variant	Disease	Class	Drug	Clinical significance	
KRAS/G12D	Colorectal Cancer	1	EXON 2 MUTATION	Cetuximab	Resistant
			MUTATION	Cetuximab	Resistant
			MUTATION	Panitumumab	Resistant
Colorectal Cancer	Colorectal Cancer	2	MUTATION	Bevacizumab	Resistant
			EXON 2 MUTATION	Cetuximab	Resistant
			G12/G13	Cetuximab	Resistant
			MUTATION	Chemotherapy	Resistant
			EXON 2 MUTATION	EGFR Inhibitor	Resistant
			EXON 2 MUTATION	Irinotecan	Sensitive
			EXON 2 MUTATION	Regorafenib	Sensitive
			MUTATION	ICM6E1000	no info
EXON 2 MUTATION	Selumetinib (AZD6244)	Sensitive			
MUTATION	Selumetinib (AZD6244)	Sensitive			
G12	(Prognostic)	no info			
Non-small Cell Lung Carcinoma	Non-small Cell Lung Carcinoma	3	MUTATION	Docetaxel	Resistant
			MUTATION	Selumetinib (AZD6244)	Resistant
Acute Leukemia	Acute Leukemia	4	G12	(Diagnostic)	Positive
Hepatocellular Carcinoma	Hepatocellular Carcinoma	4	MUTATION	Refametinib	Sensitive
			MUTATION	Sorafenib	Sensitive
Lung Adenocarcinoma	Lung Adenocarcinoma	4	MUTATION	(Prognostic)	Poor Outcome
Lung Cancer	Lung Cancer	4	G12D	Gefitinib	Resistant
			G12D	(Diagnostic)	Positive

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**CANCER 4 CAST Results**

Registration No. 20180405-99991  
Date Sample Collected 04 / 05 / 2018  
Date Sample Received 04 / 05 / 2018  
Date of Test 04 / 05 / 2018  
Date of Report 04 / 10 / 2018

**Specimen information and clinical information**  
Sample Specimen: EDTA W/B  
Medical Record:  
Family history:

**Organization & Patient Information**  
Organization: Sample  
Medical Department:  
Physician:  
Chart No.:  
Patient: Sample  
DOB (YY/MM/DD): 73/01/01  
Age / Gender: 45 / Male

**Result of Quality Control**

Specimen quality	Sequencing quality	Analysis quality
Pass	Pass	Pass

[Result]  
No pathogenic variants were detected.

**CANCER 4 CAST** Registration No. 20171111-99990

Organization	Patient	Chart No.
Sample	Sample	Sample

**Test Information**

[ 36 genes ]  
APC, ATM, BARD1, BMP1A1, BRCA1, BRCA2, BRIP1, CDH1, CDK4, CDKN2A, CHEK2, EPCAM, MEN1, MLH1, MRE11A, MSH2, MSH6, MUTHY, NBN, NFI, PALB2, PMS2, POLD1, POLE, PTEN, RAD50, RAD51C, RAD51D, RET, SDHB, SDHC, SDHD, SMAD4, STRK11, TP53, VHL  
\*TP53M: Deletion/Duplication analysis is included.

**[Limitations of the test]**  
The Cancer4Cast test aims to detect pathogenic mutations in the coding regions and adjacent regions of the genes specified above. Some exon segments (exon 9 and exons 11-15) of PMS2 with homologous genes were excluded in the analysis. No genetic rearrangements such as chromosomal aneuploidy, translocation and inversion are detected by the method of this test. Also, if there is a chromosomal aberration (mosaicism), the test result may not be accurate.  
Due to the technical limitations of the next-generation sequencing method used in this test, it is impossible to exclude the possibility of pathogenic mutations being undetected. The risk of each individual variation may depend on the type of mutation and family history. The risk of a detected mutation is a quoted value from published literature or databases based on Americans or Europeans and therefore may have ethnic differences.  
Classification and interpretation of all variations provided by this test is based on scientific information at the time of publication of this report. Classification of mutation may change if future information is updated. If there is a change in the classification of the mutation in the future, you will be notified through your doctor. Even if no pathogenic variations were identified from this test, it does not guarantee current or future health, because cancer may be caused by other factors such as age, environment and lifestyle.

**[Reference]**

- ClinVar: public archive of interpretations of clinically relevant variants. Landrum, Melissa J., et al. Nucleic acids research 44.D1 (2015): D862-D868.
- Database of Single Nucleotide Polymorphisms (dbSNP). Bethesda (MD): National Center for Biotechnology Information, National Library of Medicine. (dbSNP Build ID: 150)
- The Human Gene Mutation Database (HGMD) and its exploitation in the fields of personalized genomics and molecular evolution. Stenson, Peter D., et al. Current protocols in bioinformatics (2012): 1-13.
- Medical Genetics Information Resource (database online). Copyright: University of Washington, Seattle. 1993-2014. Available at <http://www.genetests.org> Accessed June 2, 2017.
- LabGenomics in-house variant database
- HGVS recommendations for the description of sequence variants: 2016 Update. Dunnen, Johan T., et al. Human mutation 37.6 (2016): 564-569.
- Standards and guidelines for the interpretation of sequence variants: a joint consensus recommendation of the American College of Medical Genetics (ACMG) and Genomics and the Association for Molecular Pathology, Richards, Sue, et al. Genetics in medicine (2015): 405-424.
- NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) Genetic/Familial High-Risk Assessment: Breast and Ovarian Version 2.2017. <http://www.nccn.org>
- NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) Genetic/Familial High-Risk Assessment: Colorectal Version 1.2017. <http://www.nccn.org>
- K.P. Lowy, et al. Annual screening strategies in BRCA1 and BRCA2 genetic mutation carriers: a comparative effectiveness analysis. Cancer 2012;2021-2030.

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## IMS Test Report

Inherited Metabolic disease Screening Test

Registration No.	20181227-99990
Receipt Date	2018-12-27
Analysis Date	2018-12-27
Report Date	2019-01-21

### Specimen information & Clinical opinion

Specimen: **Blood**

Specimen collection date : **2018-12-27**

Clinical opinion: **임상소견**

Family history: **가족력**

### Organization & Patient information

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### Test Result

**Pathogenic mutations associated with Glycogen storage disease Ib was detected.**

#### [Comment]

A pathogenic mutation known to be associated with Glycogen storage disease Ib was detected in the SLC37A4 gene.

Genetic counseling and confirmation tests are required for accurate diagnosis.

And a mutation known to be associated with Pompe's disease and Glycogen storage disease VII were detected in the GAA and PFKM gene.

Pompe's disease and Glycogen storage disease VII is an autosomal recessive disorder. A carrier of the recessive disorder usually have no disease-related symptoms (asymptomatic) or may appear very mild even when present. However, in the future, if the client meets a spouse with the same mutation and gives birth, there is a 25% chance of having a baby with the disease. Therefore, we recommend a genetic test for a spouse or a baby in the future. This is a screening test for genetic metabolic disorders. Genetic counseling and confirmation tests are required for accurate diagnosis.

## IMS Test Report

Inherited Metabolic disease Screening Test

Registration No. 20181227-99990

Organization	Patient	Chart No.
본사 Test	IMS	12345

Disease Type	Disease	Gene	Result	Declaration
Lysosomal storage disorders	Gaucher's disease	GBA	Not detected	Normal
	Fabry's disease	GLA	Not detected	Normal
	Pompe's disease	GAA	Not detected	Normal
	Hurler syndrome	IDUA	Not detected	Normal
	Hunter's syndrome	IDS	Not detected	Normal
	Sanfilippo syndrome A	SGSH	Not detected	Normal
	Sanfilippo syndrome B	NAGLU	Not detected	Normal
	Sanfilippo syndrome C	HGSNAT	Not detected	Normal
	Sanfilippo syndrome D	GNS	Not detected	Normal
	Morquio syndrome A	GALNS	Not detected	Normal
	Morquio syndrome B	GLB1	Not detected	Normal
	Maroteaux-Lamy syndrome	ARSB	Not detected	Normal
	Sly syndrome	GUSB	Not detected	Normal
Copper metabolic disorder	Wilson's disease	ATP7B	Not detected	Normal
Glycogen Storage Diseases	Glycogen storage disease Ia	G6PC	Not detected	Normal
	Glycogen storage disease Ib	SLC37A4	Not detected	Normal
	Glycogen storage disease III	AGL	Not detected	Normal
	Glycogen storage disease IV	GBE1	Not detected	Normal
	Glycogen storage disease V	PYGM	Not detected	Normal
	Glycogen storage disease VI	PYGL	Not detected	Normal
	Glycogen storage disease VII	PFKM	Not detected	Normal

2/3

LabGenomics Co.,Ltd.

Genetic Testing Laboratory: No. 23  
 Genetic Research Laboratory: No. 7

Medical Doctor : M.D., Ph.D.

Analysis officer : HJ Hu PH.D.

Laboratory officer : MJ Oh PH.D.

1/3

LabGenomics Co.,Ltd.

Genetic Testing Laboratory: No. 23  
 Genetic Research Laboratory: No. 7

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